

Shipping and articulator guide

Shipping

To ensure efficient and accurate processing of your order, always use the Atlantis™ CaseSafe when shipping your case to Astra Tech. Depending on your order, you will need to include specific items in the Atlantis CaseSafe.

What to include when shipping Atlantis abutment cases to Astra Tech

Cases with 1–3 adjacent edentulous spaces include:

Articulated master casts with analogs placed
Bite registration (recommended)
Model of soft tissue (recommended)

Cases with four or more edentulous spaces include:

Articulated master casts with analogs placed
Bite registration (recommended)
Model of soft tissue (recommended)
Removable diagnostic preparation that relates to the master casts

Full arch cases include:

Articulated master casts with analogs placed
Record base or occlusal bite rim (recommended)
Model of soft tissue (recommended)
Removable diagnostic preparation that relates to the master casts

Note: Case material must be disinfected prior to shipment.



Order Atlantis™ CaseSafe from your local Astra Tech representative.

Articulators

The cornerstone of any final restoration is the intraoral registration of the dental arches. A proper bite registration will ensure the best possible final result for the patient.

All cases must be articulated when shipped to Astra Tech. Make sure to verify the models are in correct occlusion. For single unit and less complicated cases an ocludator can be used.

Ceramco Vertex
Slot plastic ocludator



Ray Foster Apex or similar
Must fit in the CaseSafe™



All cases must be articulated. For larger/complicated cases the following articulators are available at Astra Tech.

Note: If using an articulator with a removable base plate listed below, it is only necessary to ship the cast mounted on the plate and not the entire articulator.

Artex
Incl. magnetic plates



Denar
Incl. screw or magnetic plates



Dentatus
Incl. screw plates



Hanau
Incl. screw plates



Kavo Protar 3
Incl. magnetic plates



Panadent
Incl. magnetic plates



Sam 3
Incl. screw or magnetic plates



Stratus 200
Incl. screw or magnetic plates



Whip Mix
Incl. screw plates

